

Visit us online now: torreypinesscientific.com

760-930-9400

Toll Free within the USA: 866-573-9104


SC20 AND SC25

MODELS SC20 DIGITAL AND SC25 PROGRAMMABLE SHAKING CHILLING/HEATING DRY BATHS

EchoTherm™ Shaking Chilling/Heating Dry Baths are ideal tools for any Molecular Biology Lab



- Shake in 2 mm orbital motion from 200 to 1000 rpm
- Peltier powered for both chilling and heating
- -10°C to 100°C (models SC20/SC25) with 50 watts of power
- -10°C to 110°C (models SC20XR/SC25XR) with 100 watts of power
- -20°C to 100°C (models SC20XT/SC25XT) with 100 watts of power
- Settable, readable, and controllable to 1°C
- Control at or near room temperature
- 5-program memory of 10-steps/program (SC25/SC25XR/SC25XT)



- Model SC20 with attachments (left)
- SC20 display (right)
- 30-day count down timer with audible alarm and user settable Auto-Off
- RS232 I/O port for computer control or data collection
- Built-in data logger
- 24 standard sample blocks available for assay plates, centrifuge tubes, test tubes, vials, and round-bottom flasks
- Custom sample blocks made on special order for most any sample container
- UL, CSA and CE compliant
- 12-month warranty
- Universal power supply provided for use anywhere in the world

Some Applications

- Enzyme reactions and deactivations
- Ligations
- Hybridizations
- PCR sample preparation
- Storing oocytes
- Storing samples at ice bucket temperatures



2713 Loker Avenue West, Carlsbad, CA 92010
Tel: 760-930-9400 | Toll Free within the USA: 866-573-9104
Fax: 760-930-9480
info@torreypinesscientific.com







760-930-9400

Toll Free within the USA: 866-573-9104

PRODUCT DATA SHEET

SC20 AND SC25

MODELS SC20 DIGITAL AND SC25 PROGRAMMABLE SHAKING CHILLING/HEATING DRY BATHS

EchoTherm™ Shaking Chilling/Heating Dry Baths are ideal tools for any Molecular Biology Lab



DIMENSIONS

7.75" (197 mm) wide x 14.5" (368 mm) front to back x 4.25" (108 mm) tall

Weight unpacked: 20 pounds (9 kgs) Weight packed: 27 pounds (12.2 kgs)