PRODUCT DATA SHEET

Visit us online now: torreypinesscientific.com

760-930-9400

Toll Free within the USA: 866-573-9104

HS55 AND HS65 ECHOTHERM[™] MODELS HS55 DIGITAL & HS65 PROGRAMMABLE 5-POSITION STIRRING HOT PLATES

()

Stir 5 vessels at the same or different digitally controlled speeds with controlled heating from room temperature to 450°C



- 12" x 12" (304 mm) ceramic glass heater top
- Control plate surface with sensor in plate or
- Control solution temperature directly with probe
- 1400 watts of heater power
- Temperature settable from room temperature to 450°C
- Heater settable and readable to 1°C
- Heater controllable to +/-1°C with PID control loop
- Temperature accuracy 1% of reading

• Temperature ramping up or down

- Platinum RTD temperature measurement system
- Model HS65 with beakers
- 5-position stirrer, all positions the same temp.
- Stir each position at the same or different speeds
- Stirrers settable and controllable to 10 rpm
- Stir from 100 to 1500 rpm

()

- 5- 800 ml beaker capacity on plate surface
- Fully programmable HS65 holds 10 programs in memory for instant recall and use
- HS65 only)
 Timer to 99 hours with audible alarm
- User settable Auto-Off
- Electronic temperature calibration traceable to NIST with certificate provided
- RS232 I/O port for control or data collection
- ETL, CSA and CE compliant
- 12-month warranty
- Available in 115VAC and 230VAC, 50/60Hz models
- HS55 simple digital control
- HS65 fully programmable control



2713 Loker Avenue West, Carlsbad, CA 92010 Tel: 760-930-9400 | Toll Free within the USA: 866-573-9104 Fax: 760-930-9480 info@torreypinesscientific.com

0/9/14 10-50 PM

Visit us online now: torreypinesscientific.com

PRODUCT DATA SHEET

760-930-9400 Toll Free within the USA:

HS55 AND HS65 ECHOTHERM[™] MODELS HS55 DIGITAL & HS65 PROGRAMMABLE 5-POSITION STIRRING HOT PLATES

۲

 \bigcirc

Stir 5 vessels at the same or different digitally controlled speeds with controlled heating from room temperature to 450°C



DIMENSIONS

۲

12.5" (318 mm) wide x 17" (432 mm) front to back x 4.5" (114 mm) tall

Weight unpacked: 24 pounds (10.9 kgs) Weight packed: 28 pounds (12.7 kgs) ۲

۲